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For *in vitro* diagnostic use

Intended use

S-Monovette® blood collection systems, needles and adapters are used together as a system for venous blood collection. S-Monovette® blood collection systems are used to collect, transport, process and store blood samples for *in vitro* diagnostic tests on serum, plasma or whole blood in the clinical laboratory. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

S-Monovette® blood collection systems consist of a plastic container, a plunger, a plunger rod and a colour-coded plastic screw cap with an integrated membrane developed especially for connecting it to an S-Monovette® (safety) needle or (Safety-) Multifly® needle, and various additives (preparations). An S-Monovette® without additive is also available. S-Monovette® blood collection systems with the additives clotting activator, clotting activator with gel, lithium heparin, lithium heparin with gel, sodium heparin, citrate, EDTA, fluoride EDTA, fluoride heparin, fluoride citrate, CPDA, as well as without additive, with additive concentrations, the volume of liquid additive and their permitted tolerances and the ratio of blood to additive complying with the requirements and recommendations of the international DIN EN ISO 6710 standard "Single-use containers for venous blood specimen collection" and the Clinical and Laboratory Standards Institute's Approved Standards (CLSI). The selection of the additive depends on the analytical test method. It is specified by the manufacturer of the test reagents and/or analysis equipment used to carry out the test. The S-Monovette® blood collection system is available as a sterile inside version or as an individually packaged version that is sterile inside and out.

Colour codes of the S-Monovette® screw caps\*:

Additive	Letter code	Cap colour DIN EN ISO 6710	Cap colour based on BS 4851 <sup>†</sup>
<b>Tube without additive</b>			
S-Monovette® neutral Z	Z	white	white
<b>Serum tube</b>			
S-Monovette® Serum CAT, with clotting activator	CAT	red	white
S-Monovette® Serum Gel CAT, with clotting activator and gel	CAT	dark yellow	dark yellow
S-Monovette® Serum Express Gel CAT, with clotting activator thrombin and gel	CAT	orange	-
<b>Heparin tube</b>			
S-Monovette® Lithium Heparin LH, with lithium heparin	LH	green	orange
S-Monovette® Lithium Heparin, liquid LH, with lithium heparin	LH	green	orange
S-Monovette® Lithium Heparin Gel LH, with lithium heparin and gel	LH	green	orange
S-Monovette® Lithium Heparin Gel* LH, with lithium heparin and rapid gel	LH	green	orange
S-Monovette® Sodium Heparin NH, with sodium heparin	NH	brown	green
S-Monovette® Ammonium Heparin AH, with ammonium heparin	n.d.	blue	blue
<b>Citrate tube</b>			
S-Monovette® Citrate 9NC, sodium citrate 3.13%	9NC	blue	green
S-Monovette® PFA, Citrate 9NC, sodium citrate 3.8%, buffered	9NC	light blue	-
<b>EDTA tube</b>			
S-Monovette® K3 EDTA K3E, with potassium EDTA	K3E	violet	red
S-Monovette® K2 EDTA K2E, with potassium EDTA	K2E	violet	red
S-Monovette® K2 EDTA Gel K2E, with potassium EDTA and gel	K2E	violet	red
<b>Tube with glycolysis inhibitors</b>			
S-Monovette® Fluoride EDTA FE, with fluoride and EDTA	FE	grey	yellow
S-Monovette® Fluoride Heparin FH, with fluoride and heparin	FH	-	yellow
S-Monovette® GlucoEXACT FC, with fluoride, citrate and EDTA	FC	pink	grey
<b>Tube for metal analysis</b>			
S-Monovette® Metal analysis LH, with lithium heparin	LH	dark blue	orange

<b>Other whole blood and plasma tubes</b>			
S-Monovette® ThromboExact, with magnesium compound	n.d.	n.d., raspberry	n.d., raspberry
S-Monovette® CTAD	n.d.	n.d., blue	n.d., green
S-Monovette® ACD	ACD	yellow	yellow
S-Monovette® CPDA	CPDA	yellow	yellow
S-Monovette® Hirudin	n.d.	n.d., olive green	n.d., olive green
S-Monovette® Homocysteine HCY-Z Gel, with gel and stabiliser	n.d.	n.d., grey	n.d., grey
S-Monovette® Homocysteine HCY-C, with citrate buffer	n.d.	n.d., violet	n.d., violet
S-Monovette® DNA Exact	n.d.	n.d., blue	n.d., blue

\* Additional screw cap colours available.

n.d.: not defined

<sup>†</sup>British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, revoked.

### SARSTEDT S-Monovette® Serum CAT / Serum Gel CAT

The S-Monovette® Serum CAT contains a synthetic granulate coated with a clotting activator (silicate) and is used for serum collection\*\*. The synthetic granulate settles between the blood clot and serum during centrifugation.

The S-Monovette® Serum Gel CAT contains a polymer-based gel, in addition to the synthetic granulate coated with a clotting activator (silicate), and is used for serum collection\*\*. The gel settles between the blood clot and serum during centrifugation and forms a stable separation layer during transport, storage and analysis.

Serum is used as sample material in routine biochemical and immunological testing and in serology.

After collecting blood, let the blood coagulate in the S-Monovette® Serum CAT / Serum Gel CAT for 30 minutes in an upright position to ensure a well-formed horizontal separation layer after centrifugation. The recommended time is based on an intact coagulation process. Blood from patients with coagulation disorders due to their illness or under anticoagulant therapy requires more time for complete coagulation.

### SARSTEDT S-Monovette® Serum Express Gel CAT

The S-Monovette® Serum Express Gel CAT contains a polymer-based gel, in addition to the synthetic granulate coated with clotting activator (thrombin), and is used for serum collection\*\*. The gel settles between the blood clot and serum during centrifugation and forms a stable separation layer during transport, storage and analysis.

After collecting blood, let the blood coagulate in the S-Monovette® Serum Express Gel for 5 minutes in an upright position to ensure a well-formed horizontal separation layer after centrifugation. The recommended time is based on an intact coagulation process, also for blood from heparinised patients. Blood from patients with coagulation disorders due to their illness or under anticoagulant therapy including high-dose heparin requires more time for complete coagulation.

Serum is used as sample material in routine biochemical and immunological testing and in serology.

### SARSTEDT S-Monovette® neutral Z

The S-Monovette® neutral Z contains no additive, which is equivalent to natural coagulation of the blood sample without clotting activators, and is used for serum collection\*\*. Centrifugation results in separation of blood clot and serum.

Serum is used as sample material for user-specific analyses.

Furthermore, the S-Monovette® neutral Z can also be used as an empty container as follows: if blood collection is to be started with a S-Monovette® Citrate in combination with a (Safety-)Multify® needle (butterfly needle), the S-Monovette® neutral Z is used as the first container (empty container) for filling the tubing of the (Safety-)Multify® needle.

### SARSTEDT S-Monovette® Lithium Heparin LH / Lithium Heparin Gel LH / Lithium Heparin Gel+ LH, Sodium Heparin NH, Ammonium Heparin AH

The S-Monovette® Lithium Heparin LH / Sodium Heparin NH / Ammonium Heparin AH contains a synthetic granulate coated with the anticoagulant lithium heparin / sodium heparin / ammonium heparin (generally 16 IU heparin/ml blood) and is used for whole blood and plasma collection\*\*. The synthetic granulate settles between the blood cells and plasma during centrifugation.

In a S-Monovette® Lithium Heparin LH without synthetic granulate, the anticoagulant lithium heparin is in a spray dose form (generally 19 IU heparin/ml blood).

The S-Monovette® Lithium Heparin Gel LH contains a polymer-based gel, in addition to the synthetic granulate coated with the anticoagulant lithium heparin (generally 25 IU heparin/ml blood), and is used for plasma collection\*\*. The gel settles between the blood cells and plasma during centrifugation and forms a stable separation layer during transport, storage and analysis.

The S-Monovette® Lithium Heparin Gel+ LH also contains a polymer-based gel that allows the centrifugation time to be reduced by approx. 50% compared with the S-Monovette® Lithium Heparin Gel.

Heparin plasma is used as sample material in routine biochemical and immunological testing. Heparin whole blood is used as sample material in routine immunological testing and in serology.

Anticoagulation with heparin occurs through activation of antithrombin.

**NOTE: No lithium measurements may be made from the S-Monovette® Lithium Heparin LH, Lithium Heparin Gel/Gel+, no sodium measurements may be made from the S-Monovette® Sodium Heparin NH and no ammonium measurements may be made from the S-Monovette® Ammonium Heparin AH. A S-Monovette® Lithium Heparin / Lithium Heparin Gel/Gel+ / Sodium Heparin / Ammonium Heparin filled with venous blood is not suitable for a blood gas analysis.**

### SARSTEDT S-Monovette® Metal analysis

The S-Monovette® Metal analysis contains the anticoagulant lithium heparin and is used for whole blood and plasma collection\*\*. Lithium heparin is in a spray dose form (19 IU heparin/ml blood).

Heparin whole blood and heparin plasma of the S-Monovette® Metal analysis is used in combination with the special safety needle 21G x 1 ½" (Art. no. 85.1162.600) as sample material for testing for trace elements/metals. The trace elements / metals were validated by atomic absorption spectroscopy (AAS). For the system safety needle and S-Monovette® Metal analysis, this results in the following maximum blank values in ng/system for these trace elements:

Element	Tl	Cd	Ni	Cr	Pb	Fe	Cu	Zn	Mn	Al	Se	Hg
ng/system	2.5	1.5	8.0	5.0	5	50	70	70	10	40	10	10

**The maximum blank values are indicated for each element on the label of the S-Monovette®.**

Trace element determinations are carried out using heparin whole blood or heparin plasma, depending on the analysis method.

Anticoagulation with heparin occurs through activation of antithrombin.

### SARSTEDT S-Monovette® Citrate 9NC

The S-Monovette® Citrate 9NC contains the anticoagulant trisodium citrate and is used for whole blood and plasma collection\*\*. Trisodium citrate is in the form of a 0.106 molar solution (3.13% trisodium citrate solution; often rounded up to 3.2%) and amounts to 10% of the nominal volume of an S-Monovette® Citrate. The mixture ratio of citrate to blood is 1:9 – 1 part by volume citrate and 9 parts by volume blood. Correct filling is absolutely necessary for the analysis.

Citrate plasma is used as sample material in routine testing in haemostaseology.

Anticoagulation occurs through calcium ion complex formation by citrate.

Please observe the NOTES under sample collection and handling for the S-Monovette® Citrate and butterfly needle combination.

### SARSTEDT S-Monovette® PFA 9NC

The S-Monovette® PFA 9NC contains the anticoagulant trisodium citrate for whole blood and plasma collection\*\*. Trisodium citrate is in the form of a 0.129 molar buffered solution (3.8% buffered; pH 5.5) and amounts to 10% of the nominal volume of a S-Monovette® PFA. The mixture ratio of citrate to blood is 1:9 – 1 part by volume citrate and 9 parts by volume blood. Correct filling is absolutely necessary for the analysis.

The whole blood of the S-Monovette® PFA 9NC is used as sample material for determination of platelet function with the SIEMENS Healthineers PFA analyser and is validated for this.

Anticoagulation occurs through calcium ion complex formation by citrate.

Please observe the NOTES under sample collection and handling for the S-Monovette® PFA and butterfly needle combination.

### SARSTEDT S-Monovette® CTAD

The S-Monovette® CTAD contains the anticoagulant trisodium citrate as well as theophylline, adenosine and dipyridamole and is used for whole blood and plasma collection\*\*. CTAD is in the form of a solution and amounts to 10% of the nominal volume of a S-Monovette® CTAD. The mixture ratio of CTAD to blood is 1:9 – 1 part by volume CTAD and 9 parts by volume blood. Correct filling is absolutely necessary for the analysis.

CTAD plasma is used as sample material for the determination of platelet factor 4 (PF4), beta(β)-thromboglobulin (βTG) and plasminogen activator inhibitor-1 (PAI-1) as well as for routine testing in haemostaseology.

Platelet activation is inhibited for a period of at least 4 hours and is used to monitor heparin therapies.

Anticoagulation occurs through calcium ion complex formation by citrate.

Please observe the NOTES under sample collection and handling for the S-Monovette® CTAD and butterfly needle combination.

**NOTE: CTAD is yellow. The S-Monovette® CTAD must be protected from artificial and natural light during storage. Exposure to strong light for more than 12 hours can deactivate the additives.**

### SARSTEDT S-Monovette® EDTA K2E / K3E

The S-Monovette® EDTA K2E / S-Monovette® EDTA K3E contains the anticoagulant K2 EDTA or K3 EDTA and is used for whole blood collection\*\*. K2 EDTA or K3 EDTA is in a spray dose form (1.6 mg EDTA/ml blood). EDTA whole blood is used as sample material for haematological analyses. Blood smears should be made within four hours after blood collection.

Anticoagulation occurs through calcium ion complex formation by EDTA.

Blood from the S-Monovette® EDTA K2E / K3E can also be used for routine immunohaematology analyses and tests for infectious diseases. The user is responsible for validating the suitability of the sample material for these analyses with the appropriate test reagents/analysis equipment incl. storage conditions.

### SARSTEDT S-Monovette® EDTA Gel K2E

The S-Monovette® EDTA Gel K2E contains a polymer-based gel, in addition to the anticoagulant K2 EDTA and is used for plasma collection\*\*. K2 EDTA is in a spray dose form (1.6 mg EDTA/ml blood). The gel settles between the blood cells and plasma during centrifugation and forms a stable separation layer during transport, storage and analysis.

EDTA plasma is used as sample material for molecular virus diagnostics.

Anticoagulation occurs through calcium ion complex formation by EDTA.

### SARSTEDT S-Monovette® Fluoride EDTA / Fluoride Heparin / Fluoride Citrate EDTA with glycolysis inhibitors

The S-Monovette® Fluoride EDTA FE contains the anticoagulant EDTA (1.2 mg/ml blood) as well as fluoride (1.0 mg/ml blood) for glycolysis inhibition and is used for plasma collection\*\*. Fluoride EDTA is in a spray dose form.

The S-Monovette® Fluoride Heparin FH contains the anticoagulant heparin (16 IU heparin/ml blood) as well as fluoride (1.0 mg/ml blood) for glycolysis inhibition and is used for plasma collection\*\*. Fluoride heparin is in a spray dose form.

Glycolysis inhibition through fluoride alone starts approx. 2 hours after filling the blood collection tube and shows its full effect after approx. 4 hours. On average, therefore, in whole blood up to 6% of glucose breaks down after a few hours and up to 10-15% after 24 hours.

The S-Monovette® GlucoEXACT FC contains the anticoagulant EDTA as well as citrate and fluoride for immediate glycolysis inhibition and is used for plasma collection\*\*. EDTA, citrate and fluoride are in the form of a solution. Because of the dilution, it is necessary to multiply by 1.16 to calculate the glucose concentration.

The S-Monovette® GlucoEXACT FC complies with the gestational diabetes guideline of the German Diabetes Association (DDG) and the national treatment guideline (NVL) for type 2 diabetes. The immediate stabilisation of the glucose concentration is maintained for up to 96 hours at room temperature in whole blood, see also SARSTEDT White Paper <https://www.sarstedt.com/en/products/diagnostic/venous-blood/s-monovetter/Literature/> "Sarstedt S-Monovette® GlucoEXACT A blood collection device for stabilizing glucose levels for 96 hours".

Fluoride EDTA (citrate) plasma is used as sample material for glucose determination.

Anticoagulation occurs through calcium ion complex formation by EDTA.

**NOTE: Fluoride can cause increased haemolysis. Further information on substances that can cause interference can be found in the respective instructions for use of the test reagent manufacturer.**

Please observe the NOTES under sample collection and handling for the S-Monovette® GlucoEXACT and butterfly needle combination.

### SARSTEDT S-Monovette® CPDA

The S-Monovette® CPDA contains the anticoagulant sodium citrate as well as citric acid, sodium phosphate, dextrose and adenine and is used for whole blood and plasma collection\*\*. CPDA is in the form of a solution. The mixture ratio of CPDA to blood is 1:6 – 1 part by volume CPDA and 6 parts by volume blood. Correct filling is absolutely necessary because of the analysis.

CPDA whole blood / plasma is used as sample material for blood group serology in transfusion medicine, for tolerability investigations (cross matching) with transfusions. Specifically, the S-Monovette® CPDA is suitable for blood group testing (ABO) including Kell and Rhesus factor, Coombs test for irregular blood group antibodies, polyspecific direct Coombs test (up to day 10), tolerability test (cross matching) with blood group compatible patient serum/plasma (major test).

Transport duration until storage is validated for a maximum of 4 hours at 2 - 21 °C. For sedimentation of the cells, the samples must first be stored in an upright position in the refrigerator for at least 16 hours. Then the aforementioned analytes can be determined for up to 35 days analogous to the usual life cycle of banked blood after storage in the refrigerator (2 - 6 °C) (with the exception of 10 days for the direct Coombs test). Suitable for the column agglutination technique e.g. on the Vision Max from Ortho Clinical Diagnostics. The sample material must not be frozen.

The suitability of the sample material for these analyses with other methods / test reagents / analysis equipment incl. storage conditions is the responsibility of the user and must be validated by the user.

Anticoagulation occurs through calcium ion complex formation by citrate.

Please observe the NOTES under sample collection and handling for the S-Monovette® CPDA and butterfly needle combination.

### SARSTEDT S-Monovette® ThromboExact

The S-Monovette® ThromboExact contains a magnesium salt as anticoagulant and is used for whole blood collection\*\*. The magnesium salt is in a spray dose form.

The whole blood of the S-Monovette® ThromboExact is used as sample material for the exclusion of false low platelet counts as a result of anticoagulant intolerance (such as EDTA, citrate, heparin), also termed pseudothrombocytopenia. Formation of platelet aggregates is prevented and determination of the actual platelet count is possible up to 12 hours after blood collection.

Inhibition of coagulation occurs through magnesium ions, with the magnesium ions inhibiting platelet-derived thromboxane A2 as well as the platelet aggregation mediated by ADP and arachidonic acid.

### SARSTEDT S-Monovette® Hirudin

The S-Monovette® Hirudin contains the anticoagulant hirudin and is used for whole blood collection\*\*. Hirudin is in a spray dose form.

Hirudin whole blood is used as sample material for the determination of platelet function on the Multiplate® apparatus (multiple platelet function analyser) from ROCHE Diagnostics and is validated for this. It is used for monitoring therapy with platelet inhibiting medication and for the proof or exclusion of platelet function disorders.

The S-Monovette® Hirudin was developed especially for determining platelet function in cooperation with the company Verum Diagnostica, now ROCHE Diagnostics, on the Multiplate® apparatus.

Inhibition of coagulation occurs through hirudin by direct inhibition of thrombin, thus allowing platelet function diagnostics in the native condition.

### SARSTEDT S-Monovette® Homocysteine HCY-Z Gel / HCY-C

The S-Monovette® Homocysteine HCY-Z Gel contains a special stabiliser and a polymer-based gel, in addition to the synthetic granulate coated with a clotting activator (silicate), and is used for serum collection\*\*. The gel settles between the blood clot and serum during centrifugation and forms a stable separation layer during transport, storage and analysis. The additive is in a spray dose form.

Homocysteine is stabilised for up to 8 hours in whole blood at room temperature. If centrifugation takes place within the first 8 hours with a correctly formed separation layer between the blood clot and serum, homocysteine is stabilised for up to 96 hours.

The serum of the S-Monovette® Homocysteine HCY-Z-Gel is used as sample material for homocysteine determination. The S-Monovette® Homocysteine HCY-C contains acid citrate (pH 4.3) as anticoagulant and is used for plasma collection\*\*. Citrate is in the form of a solution. Because of the dilution, it is necessary to multiply by 1.11 to calculate the homocysteine concentration.

Homocysteine in whole blood is stabilised for up to 6 hours at room temperature and up to 48 hours at 4°C.

The plasma of the S-Monovette® Homocysteine HCY-C is used as sample material for homocysteine determination.

Inhibition of coagulation in the S-Monovette® Homocysteine HCY-C occurs through calcium ion complex formation by citrate.

Please observe the NOTES under sample collection and handling for the S-Monovette® HCY-C and butterfly needle combination.

### SARSTEDT S-Monovette® DNA Exact

The S-Monovette® DNA Exact contains an anticoagulant based on EDTA and is used for whole blood collection\*\*. EDTA is in a spray dose form.

EDTA whole blood is used as sample material for stabilisation of gDNA. The stabilisation property for gDNA is validated for 5 days at 35 °C, 14 days at room temperature (22 °C), 28 days under refrigeration (4 °C) and at least 1 year at -20 °C. In addition, the stabilised specimens can be frozen and thawed for 5 cycles without loss of quality, see also SARSTEDT White Paper <https://www.sarstedt.com/en/products/diagnostic/venous-blood/s-monovetter/Literature/> "SARSTEDT S-Monovette® DNA Exact – A New IVD Certified Blood Collection Tube for Collection, Transport and Stabilization of Whole Blood for Genomic DNA Analysis".

Anticoagulation occurs through calcium ion complex formation by EDTA.

\*\*Suitability of the sample material depends on the analyte and test reagent / analysis equipment. Corresponding recommendations of the manufacturer must be observed.

## Safety information / Warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from blood and possible exposure to potentially infectious specimen material and transmissible pathogens.
2. Treat all biological specimens and sharp blood collection utensils (needles, adapters) according to the guidelines and procedures in your institution. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Use safety needles / Safety-Multifly® needles with integrated needle guard. Always follow the safety guidelines and procedures of your facility.
3. Applies to individually sterile-packed products: Check each individual package for damage before use and do not use if packaging is damaged.
4. The product is intended for single use. Dispose of all sharp objects (needles, adapters) for blood collection in suitable containers for hazardous biological waste.
5. If blood is collected from an intravenous (IV) line, it must be ensured that the line is properly flushed in accordance with the institution (=cleaned of the IV solution) before beginning blood collection. Properly flushing the line prevents faulty analysis results.
6. Underfilling the S-Monovette® leads to an incorrect ratio of blood to preparation/additive and can lead to false analysis results.
7. The blood collected with the S-Monovette® and processed is not intended to be re-injected into the human body.
8. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.

## Storage

Store the product at room temperature until use.

## Limitations

1. The storage period and temperature of a filled S-Monovette® depend on the stability of the analytes for investigation. The assessment is made by the laboratory, or information is obtained from the manufacturer's instructions for use for the test reagents / analysis device.
2. Even if plasma or serum is separated by centrifuging a S-Monovette® and/or a barrier is present, all cells are not necessarily completely separated. Residual metabolism or natural degradation can affect the concentrations of analytes. The stability of the analytes should be evaluated with respect to the storage container and the conditions of the respective laboratory.
3. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use for the test reagents / analysis device.
4. Analyses for trace elements / metals are to be conducted with the S-Monovette® Metal Analysis + Safety Needle Metal Analysis system.

## Specimen collection and handling

**READ ALL OF THIS DOCUMENT BEFORE STARTING TO DRAW VENOUS BLOOD.**

### Equipment required for specimen collection

1. All necessary S-Monovette® blood collection systems.
2. (Safety) needles or (Safety-)Multifly® needles
3. Gloves, gown, goggles or other suitable protective gear to protect from bloodborne pathogens or potentially infectious materials.
4. Labels for patient identification.
5. Disinfectant to clean the collection site (follow guidelines of the facility for sterile specimen collection to prepare the collection site). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean, disposable swabs.
7. A tourniquet.
8. Plaster or dressing.
9. Sharps disposal container for the safe disposal of used material.

### Recommended order of collection:

1. Blood culture vials
2. Tubes with citrate or clotting activator (with/without gel, Serum CAT / Serum Gel CAT)
3. Tubes with clotting activator or citrate
4. Tubes with heparin with/without gel
5. Tubes with EDTA with/without gel
6. Tubes with glycolysis inhibitors
7. Tubes with other additives

**NOTE:** *When starting collection with a citrate tube or another tube with a liquid additive such as PFA, CTAD, CPDA, HCY-C or GlucoEXACT in combination with a butterfly needle (e.g. (Safety-)Multifly® needle), collection with an empty tube (e.g. S-Monovette® neutral Z) first is recommended to avoid underfilling of the tube due to filling of the tubing of the butterfly needle. This can ensure the correct filling of the tube and thus a correct mixing ratio (liquid to blood).*

**NOTE:** *Follow the regulations of your institution for the order of collection.*

### Preventing reflux

Most S-Monovettes contain chemical additives. To prevent reflux, the following should be noted:

1. Use only the materials recommended in the instructions for use.
2. Position the patient so that there is easy access to the vein, and the arm or other puncture site on the patient is directed downwards.
3. Hold the S-Monovette® with the screw cap pointing up and when using a (Safety-)Multifly® needle, also below the level of the arm / puncture site.
4. It must be ensured that the contents of the S-Monovette® do not touch the screw cap or the end of the needle during the venous puncture.
5. Loosen the tourniquet while the blood is flowing into the S-Monovette®.

## Freezing/thawing

According to the recommendations of the WHO (WHO/DIL/Lab/99.1 Rev.02), it is recommended to separate serum/plasma from blood cells before freezing. The following recommendations for freezing apply to the S-Monovette®:

### Freezing below 0 °C

- Use at temperatures below freezing must be tested under real conditions, as the stress on the containers can vary widely depending on the medium to be frozen and the routine conditions (see section on Material characteristics in the SARSTEDT catalogue).
- The strength values of plastics are reduced in the temperature range below 0 °C. Mechanical stress should therefore be avoided.
- The conditions for freezing must be selected so that the contents of the S-Monovette® freeze evenly or from bottom to top. The S-Monovette® should have enough room in a holder or storage box to allow it to expand. Holders made of polystyrene or metal are unsuitable as they can lead to expansion tears.

### Freezing to -20 °C – Check before freezing:

- Interferences: Check whether freezing will interfere with the blood specimen or analysis (e.g. haemolysis, analyte stability).
- Centrifugation: Necessary depending on the planned analysis. Take note of centrifugation conditions, see this document under "Centrifugation" / SARSTEDT website <https://www.sarstedt.com/en/service/centrifugation/>, and perform a visual specimen control.
- Haemolysis: Whole blood is haemolysed when freezing: To prevent this, either the S-Monovette® Gel, or after centrifugation, the Seraplas® valve filter should be used or the specimen should be pipetted into a secondary container.

### Freezing to -20 °C – Implementation:

- Cooling down: Cool down the S-Monovette® while it is standing upright for a period of 45 – 60 min. from room temperature to +4 °C before it can be frozen at -20 °C.
- Thawing: Allow to thaw at least 45 min. at room temperature while standing upright. Mechanical stress must also be prevented in this step. Thawing too quickly can impair the analysis results.

### Freezing and deep-freezing below -20 °C

- Deep-freezing to below -20 °C has not been tested by the manufacturer. Due to the number of possible influencing factors, it is recommended to conduct freezing tests under routine laboratory conditions.

### Freezing the S-Monovette® Gel:

- Freezing is implemented as under the section – Freezing to -20 °C – Implementation.
- It is known that the gel layer can be changed due to the freezing process.
- Due to the number of possible influencing factors, it is recommended to conduct freezing tests under routine laboratory conditions.
- **Thawing:** Allow the specimen to thaw while standing upright. Transfer the specimen material from the primary to a secondary container after thawing by pipetting (not decanting) it up to approx. 2 mm above the gel layer. Discard the rest.

**NOTE:** For the stability of the analytes, follow the instructions for use of the test reagents / analysis device manufacturer.

## Blood collection

### Collection techniques

A distinction is made between two collection techniques: the aspiration technique and the vacuum technique.

The aspiration technique enables gentle blood collection with a continuous, slow blood flow, by controlled pulling on the plunger rod. In this way, the blood flow can be adapted immediately to all vein conditions and circumstances as well as reducing haemolysis.

The vacuum technique using a pre-evacuated tube allows blood to be collected according to the negative pressure principle with a continuous rapid blood flow. This technique is thus ideal for good vein status and easy collection conditions.

Handling of blood collection: see also the handling video for the S-Monovette® in the aspiration technique or vacuum technique with the safety needle or Safety-Multifly® needle:

[www.sarstedt.com/en/download/videos/Diagnostics/Venous\\_blood](http://www.sarstedt.com/en/download/videos/Diagnostics/Venous_blood)

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING BLOOD COLLECTION.


1. Select a suitable S-Monovette® for the required specimen.
2. Select a site for the venipuncture. Please follow the guidelines of the institution for preparing the venipuncture site.

### Blood collection with the aspiration technique

**Important! Use S-Monovette® only for blood collection, not for injection.**

Use only with these needles (S-Monovette® (safety) needles, or (Safety-)Multifly® needles/adapters for the S-Monovette®).

Use with a safety needle is shown here:

1. Immediately prior to the venipuncture, put the S-Monovette® into the integrated holder of the safety needle **A** and secure it by **twisting clockwise slightly B**. Apply a venous tourniquet if necessary (max. 1 minute). Prepare the puncture site with a suitable disinfectant. Do not touch the puncture site after disinfection. Position the patient so that there is easy access to the vein, and the arm or other puncture site on the patient is directed downwards.

2. Remove the needle guard. Puncture the vein, slowly retract the plunger and loosen the tourniquet while the blood is flowing into the S-Monovette® **C**. Wait until the flow of blood stops to allow correct filling.

3. Move the S-Monovette® by **twisting it anticlockwise slightly D** and then pull (disconnect) it from the safety needle **E**. The safety needle remains in the vein. For multiple specimens – as described above – the subsequent S-Monovettes are connected to the safety needle and additional specimens are taken. Completing blood collection: First disconnect the S-Monovette® **D + E**, then withdraw the safety needle from the vein. Using one hand, either activate the needle guard on a stable, level surface – until the needle snaps into the needle guard with a **perceptible and audible "click"** – or by pressing the index finger against the lower end of the needle guard, see also instructions for use for the S-Monovette® (safety) needle or (Safety-)Multifly® needle. Invert each S-Monovette® once immediately after disconnecting it from the safety needle, and after completing blood collection, carefully invert all specimens with anticoagulants several times. For transport and centrifugation, retract the plunger rod until the plunger snaps audibly into place **F** (plunger clicks in the bottom of the S-Monovette®) and break off the plunger rod **G**.

1. Press a clean, dry swab against the puncture site until bleeding stops.
2. As soon as coagulation has set in, you may apply a dressing if desired.
3. Dispose of the safety needle in a suitable container for hazardous biological waste.


## Blood collection with the vacuum technique

**Important! Use S-Monovette® only for blood collection, not for injection.**

Use only with these needles (S-Monovette® (safety) needles, or (Safety-)Multifly® needles/adapters for the S-Monovette®).

Use with a safety needle is shown here:

- 1** To begin blood collection gently, SARSTEDT recommends collecting the first S-Monovette® using the aspiration technique, see VENIPUNCTURE TECHNIQUE Aspiration technique **A – E**. Then continue with the vacuum technique.

Before collecting blood, use the plunger rod to pull the plunger to the bottom of the S-Monovette® and allow it to snap into place audibly. **A** Then the plunger rod is broken off **E**.

- 2** Slide this prepared S-Monovette® into the integrated holder of the S-Monovette® safety needle and secure it by **twisting it clockwise slightly C + D**.

- 3** Wait until the flow of blood stops to allow correct filling.

- 4** Move the S-Monovette® by **twisting it anticlockwise slightly E** and then pull (disconnect) it from the safety needle **F**. The safety needle remains in the vein.

For multiple specimens – as described above – the subsequent S-Monovettes are connected to the safety needle and additional specimens are taken.

Completing blood collection:

First disconnect the S-Monovette® **E + F**, then withdraw the safety needle from the vein. Using one hand, activate the needle guard either on a stable, level surface – until the needle snaps into the needle guard with a **perceptible and audible "click"** – or by pressing the index finger against the lower end of the needle guard.

Invert each S-Monovette® once immediately after disconnecting it from the safety needle, and after completing blood collection, carefully invert all specimens with anticoagulants several times.

1. Press a clean, dry swab against the puncture site until bleeding stops.
2. As soon as coagulation has set in, you may apply a dressing if desired.
3. Dispose of the safety needle in a suitable container for hazardous biological waste.

## Centrifugation

IMPORTANT! SARSTEDT S-Monovettes are designed for 4,000 x g. S-Monovettes with a diameter of 8 mm (S-Monovette® Paediatrics), are an exception, and are currently validated up to 2,500 x g. Only use suitable tube holders or inserts. Centrifuging S-Monovettes with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the S-Monovettes, which can release potentially infectious material.

Centrifuge inserts are to be selected according to the size of the S-Monovettes used. The relative centrifugal force has the following ratio to the selected rpm:

$RCF = 11.2 \times r \times (\text{rpm}/1,000)^2$ ,

RCF: "Relative centrifugal force", also called "g-force",

"rpm": "Revolutions per minute",

"r": "Radius of rotation from the middle of the centrifuge to the bottom" of the S-Monovette®, (cm).

S-Monovettes without gel can be centrifuged in centrifuges with a fixed-angle or swing-out rotor.

S-Monovettes with gel are designed solely for use in centrifuges with swing-out rotor. Centrifugation in centrifuges with a fixed-angle rotor is not validated by SARSTEDT and is not recommended.

S-Monovettes should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the users themselves.

It must be ensured that the S-Monovettes fit well in the centrifuge inserts. S-Monovettes that protrude above the insert can get caught in the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

**CAUTION! Do not remove broken S-Monovettes by hand.**

Refer to the centrifuge instructions for use for information on disinfecting the centrifuge.

Based on BS 4851* (EU code)	DIN EN ISO 6710	S-Monovette®	Relative centrifugal force (g)				
			2000 x g	2500 x g	3000 x g*	3500 x g*	4000 x g*
		Serum	10 min	10 min	6 min	4 min	4 min
		Serum Gel	15 min	10 min	4 min	4 min	4 min
—		Serum Express Gel	15 min	10 min	4 min	n.v.	n.v.
		Lithium Heparin	10 min	10 min	7 min	7 min	7 min
		Lithium Heparin Gel	15 min	15 min	10 min	7 min	7 min
		Lithium Heparin Gel+	8 min	7 min	5 min	4 min	4 min
		EDTA	n.v.	n.v.	7 min	6 min	5 min
		EDTA Gel	15 min	10 min	10 min	7 min	7 min
		Citrate	9 min	8 min	7 min	6 min	5 min
		EDTA, Fluoride Heparin	9 min	8 min	7 min	6 min	5 min
		GlucoEXACT	9 min	8 min	7 min	6 min	5 min
		Citrat PBM 1.8 ml Centrifuge radius > 17 cm	9 min	8 min	7 min	6 min	5 min
		Citrate PBM 1.8 ml Centrifuge radius > 9 to ≤ 17 cm	n.v.	n.v.	10 min	n.v.	n.v.

n.v. = not validated

\* Applies to all S-Monovettes except with 8 mm diameter (S-Monovettes Paediatrics)

Centrifugation at 20°C

\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, revoked.

## Information on separation material

The flow properties of the separation material (gel) are temperature dependent. If S-Monovettes are cooled before or during centrifugation, a continuous layer of the separation gel is no longer ensured.

S-Monovettes must not be centrifuged again / re-centrifuged.

## Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

## Product-specific standards and guidelines as amended/updated

CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.

CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

CLSI\* H21 "Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

### Standard for containers for venous blood collection:

DIN EN ISO 6710: Single-use containers for human venous blood specimen collection

\*CLSI (Clinical and Laboratory Standards Institute)

## Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture
	Sterile fluid path
	Do not re-sterilise

### Additionally applies to individually sterile-packed products:

	Do not use if the packaging is damaged
	Single sterile barrier system with external protective packaging
	Sterilisation by irradiation

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, отменен.

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CLSI\* GP39 „Tubes and Additives for Venous and Capillary Blood Specimen Collection“, Approved Standard.  
 CLSI\* GP41 „Collection of Diagnostic Venous Blood Specimens“, Approved Standard.  
 CLSI\* GP44 „Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests“, Approved Guideline.  
 CLSI\* H21 „Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays“, Approved Guideline.  
 WHO/DIL/LAB/99.1 Rev02 „WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002“.

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- CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
- CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.
- CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
- CLSI\* H21 "Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline.
- WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

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\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, wycofano.

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CLSI\* GP39 „Tubes and Additives for Venous and Capillary Blood Specimen Collection”, Approved Standard.

CLSI\* GP41 „Collection of Diagnostic Venous Blood Specimens”, Approved Standard.

CLSI\* GP44 „Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests”, Approved Guideline.

CLSI\* H21 „Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays”, Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 „WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002”.

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CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.  
 CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.  
 CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.  
 CLSI\* H21 "Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline.  
 WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

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- CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard
- CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard
- CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline
- CLSI\* H21 "Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline
- WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

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CLSI\* GP39 " Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI\* GP41 " Collection of Diagnostic Venous Blood Specimens", Approved Standard.

CLSI\* GP44 " Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

CLSI\* H21 " Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 " WORLD HEALTH ORGANIZATION, et al.Use of anticoagulants in diagnostic laboratory investigations.Geneva:World Health Organization, 2002".

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