	WORK INSTRUCTION FM-DVSOP2002-01 Rev. 001 Eff. Date: 25-OCT-2016 CO No.: CO-10029546		Doc. Number AS10016053-002	Revision 100
	WI Type: ARTWORK SPECIFICATION		Effective Date: 11-APR-2019	Page 1 of 33
Title: INSTRUCTIONS FOR USE – BLUE LINE SILICONISED TRACHEOSTOMY TUBES (10016053-002)				

DOCUMENT VERIFICATION

Prior to using this document, the user is responsible for verifying that the revision and effective date are current.

REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MTIJ		
Rev.	Effective Date	Change Order No.	Changes Made to Document	
100	11-APR-2019	CO-10077589	PWR: 5508 This artwork replaces 10016053-001 (AS10016053-001 Revision 100). Updates include the following changes; <ul style="list-style-type: none"> • Legal Manufacturer address changed from Hythe, UK to Nathan Lane, Minneapolis. • EU Representative address changed from Ashford, UK to Czech Republic address. • New Notified Body, changed from CE0086 to CE2797. • Copyright date updated to 2019 in Trademark Statement. 	

LANGUAGES

en	fr	de	it	es	pt	nl	sv	da	no
English	French	German	Italian	Spanish	Portuguese	Dutch	Swedish	Danish	Norwegian
Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
fi	el	ja	cs	sk	hu	pl	sl	et	lv
Finnish	Greek	Japanese	Czech	Slovakian	Hungarian	Polish	Slovenian	Estonian	Latvian
Y	Y	N	Y	N	Y	Y	N	N	N
lt	tr	bg	ro	ru	hr	zh	ko	ar	
Lithuanian	Turkish	Bulgarian	Romanian	Russian	Croatian	Chinese	Korean	Arabic	
N	Y	Y	Y	N	N	N	N	N	

ARTWORK DESIGN NOTES

Part Number: 10016053-002




Size: 150mm (W) x 210mm (H). **Finished Folded Size:** 75mm (W) x 210mm (H).

IFU Format: Printed Instruction Booklet, 60 pages.

- A) For detailed IFU specifications refer to the purchase specification.
- B) Artwork Content and Specification Approval via electronic approval.
- C) Approval of this document is for final format and layout only. All translations have been previously approved as part of the approval and release process of AS10016053-001 Revision 100.

	WORK INSTRUCTION FM-DVSOP2002-01 Rev. 001 Eff. Date: 25-OCT-2016 CO No.: CO-10029546		Doc. Number AS10016053-002	Revision 100
	WI Type: ARTWORK SPECIFICATION		Effective Date: 11-APR-2019	Page 2 of 33
Title: INSTRUCTIONS FOR USE – BLUE LINE SILICONISED TRACHEOSTOMY TUBES (10016053-002)				

FOLDING DIAGRAM

	WORK INSTRUCTION FM-DVSOP2002-01 Rev. 001 Eff. Date: 25-OCT-2016 CO No.: CO-10029546	Doc. Number AS10016053-002	Revision 100
		WI Type: ARTWORK SPECIFICATION	Effective Date: 11-APR-2019
Title: INSTRUCTIONS FOR USE – BLUE LINE SILICONISED TRACHEOSTOMY TUBES (10016053-002)			

IFU ARTWORK

(To scale)

smiths medical
portex[®]

Blue Line[®] Siliconised PVC Tracheostomy Tube

- en** Instructions for use
- fr** Mode d'emploi
- de** Gebrauchsanweisung
- it** Istruzioni per l'uso
- es** Instrucciones de uso
- pt** Instruções de utilização
- nl** Gebruiksaanwijzing
- sv** Bruksanvisning
- da** Brugsanvisning
- no** Bruksanvisning
- fi** Käyttöohjeet
- el** Οδηγίες χρήσης
- cs** Návod k použití
- hu** Használati utasítás
- pl** Instrukcja stosowania
- tr** Kullanım Talimatı
- bg** Указания за употреба
- ro** Instrucțiuni de utilizare

Manufacturer:
Smiths Medical ASD, Inc.
 6000 Nathan Lane North
 Minneapolis, MN 55442 USA
 Tel: 1 800 258 5361 (US/CA)
 Tel: +1 614 210 7300

European Representative:
Smiths Medical Czech Republic a.s.
 Olomoucká 306, Hranice 1 - Město,
 753 01 Hranice, Czech Republic
 Tel: +44 (0) 1233 722100

www.smiths-medical.com

2797

10016053-002_2019-03

English

Portex® Blue Line® Siliconised PVC Tracheostomy Tube

INSTRUCTIONS FOR USE

These instructions apply to the use of the following Portex® products:

100/506/030-090	Uncuffed with 15mm connector
100/510/060-100	Traditional cuffed with 15mm connector
100/518/060-100	Profile™ Low Pressure/High Volume cuffed with 15mm connector
100/535/030-100	Uncuffed, single fenestration without 15mm connector
100/536/030-100	Uncuffed, single fenestration with 15mm connector
100/537/060-100	Uncuffed, double fenestration with 15mm connector
109/5xx/xxx	Products with the prefix 109 that are manufactured to order

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

DESCRIPTION

A sterile range of single use Tracheostomy Tubes manufactured from Siliconised Polyvinyl Chloride (PVC) incorporating the following features:

- Thermosensitive material with sufficient initial rigidity for intubation which then conforms to the individual patient's upper respiratory tract at body temperature ensuring minimum trauma.
- Anatomically shaped design with radio-opaque Blue Line® to confirm correct tube placement by X-Ray.
- Cuffed tracheostomy tubes are intended for patients at risk of aspiration and/or requiring manual or mechanical ventilation.
- Uncuffed tracheostomy tubes are intended for paediatric and adult spontaneously breathing patients, allowing speech and access to the trachea for suctioning. The 15mm connector allows attachment of humidification devices (Thermovent® T 100/570/015). Uncuffed tracheostomy tubes without a 15mm connector are intended for patients with a permanent tracheostomy. The low profile flange is ideal for wearing under clothing and tracheostomy bibs.
- Uncuffed fenestrated tracheostomy tubes allow, in addition to the above, the patient to breathe through the fenestrations as well as around the outside of the tube reducing the effort and work of breathing necessary.

Cuff Resting Diameter (mm)

Tracheostomy Tube I.D.	6.0	7.0	7.5	8.0	9.0	10.0
Traditional Cuff	11	14	14	16	18	20
Profile™ Cuff	19	22	24	27	30	32

INDICATIONS FOR USE

Portex® Blue Line® Siliconised Tracheostomy Tubes are indicated for airway maintenance of tracheostomised patients.

CONTRAINDICATIONS

Never use uncuffed fenestrated tubes in patients at risk of aspiration or who may require manual or mechanical ventilation, because material may enter the lungs and ventilation gases will leak through both the fenestrations and around the outside of the tracheostomy tube.

ADVERSE REACTIONS

Reported adverse reactions associated with tracheostomy tubes are many and diverse. Consult standard textbooks and medical literature for specific adverse reaction information.

WARNINGS

- Prior to removal of cuffed Tracheostomy Tubes, all air must be completely removed from the cuff to prevent damage to trachea and stoma.
- If the Tracheostomy Tube is lubricated prior to insertion, ensure the lubricant does not occlude the lumen of the tube preventing ventilation of the patient.
- Do not inflate with a measured volume of air, or by feel of pressure from the syringe, since little resistance should be felt during inflation.
- During anaesthesia, nitrous oxide may diffuse into the cuff causing an increase or decrease in cuff pressure.
- Cuff pressure should be adjusted and monitored routinely. Over-inflation of the cuff may result in permanent damage to the trachea.
- Contact with electrosurgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air or ignite in an enriched oxygen environment (e.g. anaesthesia).
- Never use uncuffed fenestrated tubes in patients at risk of aspiration or who may require manual or mechanical ventilation, because material may enter the lungs through the fenestrations and ventilation gases will leak through both the fenestrations and around the outside of the tracheostomy tube.
- The position and patency of the fenestration(s) in the Portex® Fenestrated Tracheostomy Tube should be checked immediately after insertion to determine optimal position of the fenestra, and frequently thereafter to ensure that tissue granulation has not caused an obstruction.
- Do not use this tracheostomy tube in patients with abnormal upper airway anatomy or pathology as this may result in partial or total airway obstruction.
- Do not use suction catheters with single fenestration tracheostomy tubes to prevent the catheter from exiting the fenestra and possibly damaging the tracheal mucosa or tracheostomy tube cuff.
- During and after attachment of the breathing system to the tracheostomy tube connector, avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion.
- Tracheostomy tubes must be changed regularly to suit individual patient's needs.
- Do not use a Portex® Purple Decannulation Cap with the above products because it may be difficult to remove in a critical or emergency situation, potentially resulting in obstruction of the patient's airway.
- If a serial percutaneous tracheostomy dilator is used as a loading dilator/obturator to insert a cuffed tracheostomy tube, ensure the dilator can be inserted and removed from the tracheostomy tube without undue force before use.
- This product contains DEHP – data derived from animal and epidemiological studies outline the following risks, but are not limited to: Adverse effects on fetal development, kidney and liver toxicity, negative impact on fertility. Special attention on the application of phthalate products should be given for the following patient populations: paediatrics, pregnant and nursing women. For further information on phthalates please refer to the Smiths Medical website: www.smiths-medical.com/phthalates

PRECAUTIONS

- The security of all breathing system connectors should be checked when the circuit is established and frequently thereafter. Disconnection may be facilitated with the use of a disconnection wedge (100/555/100).
- Patients should be adequately humidified to minimise encrustation of the tracheostomy tube lumen and prevent mucosal damage.
- The patency of the Tracheostomy Tube lumen must be assured by regular suctioning. Check routinely and replace as required to maintain a patent airway. Maximum recommended period of use 30 days.
- Devices used in or during inflation of the cuff must be clean and free from all foreign matter. The inflation device should be removed from the inflation valve immediately after use, and the dust cap fitted.
- Guard against cuff damage by avoiding contact with sharp edges.
- The inflation line valve may interfere with Magnetic Resonance Imaging (MRI) picture clarity. Ensure the valve is positioned away from area being scanned.
- If Portex® tracheostomy tubes are used outside the hospital, the patient must be instructed by a health care professional in the safe use and handling of the product. A patient care booklet is also available from Smiths Medical Customer Services.

- 8 Only sterile saline should be used to clean accessible parts of the tube whilst *in situ*. If the tracheostomy tube is removed from the patient, it must be discarded.
- 9 Blue Line® Siliconised Tracheostomy Tubes are designed and intended for single use. They should not be re-sterilised by any method.
- 10 Repositioning of the *in situ* tracheostomy tube while the cuff is inflated should be avoided.

INSTRUCTIONS FOR USE

- 1 The integrity of the cuff and inflation system should be checked prior to insertion.
- 2 Prior to use, unclip the flange from the obturator (where fitted) and check that the obturator can be removed from the tube. Re-insert the obturator.
- 3 A small amount of water soluble lubricant can be applied to the tube and obturator tip to ease insertion.
- 4 Insert tracheostomy tube in accordance with currently accepted medical techniques.
- 5 Inflate the cuff (where fitted) with the minimum amount of air to provide an effective seal.
Note: Should the product cuff remain partially stuck to the tube during inflation, gentle manipulation will help to release the cuff.
- 6 Secure the tracheostomy tube with the cotton tapes provided or a Portex® Tracheostomy Tube Holder (100/503/200).
- 7 If a Portex® Fenestrated Tracheostomy Tube is used, attach the supplied label in a prominent place to alert clinical staff to the presence of a fenestrated tracheostomy tube.
- 8 If a patient requires manual or mechanical ventilation whilst using a Fenestrated Tracheostomy Tube, a suitable tube and non-fenestrated inner cannula must be in place. We recommend the tube be replaced with a Portex® Blue Line Ultra® Fenestrated Tracheostomy Tube using a Portex® Tracheostomy Change Kit (Ref: 100/822 or 100/823).

Profile, Thermovent, Blue Line, Blue Line Ultra, Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2019 Smiths Medical. All rights reserved.

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en English	fr Français	de Deutsch	it Italiano	es Español	pt Português
 Caution Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 0434A	Attention	 Vorsicht	Attenzione	Precaución	Precaução
 Do Not Reuse Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. ISO – 15221 Symbols – Used with Medical Device Labelling. ISO 7000: 1051	Ne pas réutiliser	 Nicht zur Wiederverwendung	Non riutilizzare	No reutilizar	Não reutilizar
 Catalogue Number Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 2493	Numéro de catalogue	 Bestellnummer	Numero di catalogo	Número de catálogo	Número de catálogo
 Batch Code Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 2492	Numéro de lot	 Chargenbezeichnung	Codice del lotto	Código de lote	Código de Lote
 Date of Manufacture Indicates the date when the medical device was manufactured YYYY-MM-DD. ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 2497	Date de fabrication	 Herstellungsdatum	Data di fabbricazione	Fecha de fabricación	Data de fabrico
 Use by Indicates the date after which the medical device is not to be used YYYY-MM-DD. ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 2607	Utiliser avant le	 Verwendbar bis	Utilizzare entro	Fecha de caducidad	Utilizar até
 Manufacturer Indicates the medical device manufacturer. ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 3082	Fabricant	 Hersteller	Fabbricante	Fabricante	Fabricante
 Not made with natural rubber latex Large X indicates natural rubber latex (NRL) is not used as a material in the medical product or product container. ISO-15223-1 Symbols-Used with Medical Device Labelling Annex B.2. ISO 7000: None	Fabriqué sans latex naturel (caoutchouc)	 Nicht mit Latex aus Naturkautschuk hergestellt	Prodotto senza lattice di gomma naturale	Fabricado sin látex de caucho natural	Não fabricado com látex de borracha natural
 Contains or presence of phthalate Contains or presence of phthalate: bis (2-ethylhexyl) phthalate (DEHP). EN 15986:2011 Symbols – Used with Medical Device Labelling – Requirements for labelling of medical devices containing phthalates. Derived from ISO 7000: 2725 www.smiths-medical.com/phthalates	Contient du phtalate ou présence de phtalate : phtalate de bis(2-éthylhexyle) (DEHP)	 Enthält oder Anwesenheit von Phthalaten: Diethylhexylphthalat (DEHP)	Contiene o presenza di ftalati: di(2-etilhexil) ftalato (DEHP)	Contenido o presencia de ftalato: di(2-etilhexil) ftalato (DEHP)	Conteúdo ou presença de ftalato: ftalato de bis(2-etilhexilo) (DEHP)
 Package Quantity Indicates number of pieces in package. Symbol is read and understood by ordinary individuals. ISO 7000: None	Quantité	 Menge	Quantità	Cantidad	Quantidade
 Cuff Resting Diameter Indicates Resting Diameter of Cuff as determined by EN 1282-1 & 2:1997 Annex B, where XX is the numerical value. Symbol is read and understood by ordinary individuals. ISO-7000: None	Diamètre du ballonnet au repos	 Cuffdurchmesser (nicht gefüllt)	Diametro cuffia sgonfia	Diámetro del balón en reposo	Diâmetro da braçadeira em repouso
 Authorized Representative in the European Community Indicates the Authorized representative in the European Community. ISO-15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: None	Mandataire dans la Communauté européenne	 Autorisierter Vertreter in der Europäischen Gemeinschaft	Mandatario per la Comunità Europea	Representante autorizado en la Comunidad Europea	Representante Autorizado na Comunidade Europeia
 Sterilized using ethylene oxide Indicates a medical device that has been sterilized using ethylene oxide ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 2501	Stérilisé à l'oxyde d'éthylène	 Sterilisation mit Ethylenoxid	Metodo di sterilizzazione con ossido di etilene	Método de esterilización utilizando óxido de etileno	Esterilizado utilizando óxido de etileno
 Do not use if packaged is damaged Indicates a medical device that should not be used if the package has been damaged or opened ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 2606	Ne pas utiliser si l'emballage est endommagé.	 Inhalt bei beschädigter Verpackung nicht verwenden.	Non utilizzare se la confezione è danneggiata.	No utilizar si el envase está dañado.	Não utilizar se a embalagem estiver danificada.
 Keep dry Indicates a medical device that needs to be protected from moisture. ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 0626	Conservar au sec	 Trocken aufbewahren	Tenere all'asciutto	Manténgase en lugar seco	Manter seco
 Keep away from sunlight Indicates a medical device that needs protection from light sources ISO – 15223-1 Symbols – Used with Medical Device Labelling. ISO 7000: 0624	Conservar à l'abri de la lumière du soleil	 Vor Sonnenlicht schützen	Tenere al riparo dalla luce solare	Manténgase lejos de la luz solar	Manter afastado da luz solar

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Profile, Thermovent, Blue Line, Blue Line Ultra, Portex, designmerket Portex og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet ® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.