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Rectal catheter

Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PVC
Compatibility:	With lavage syringe Janette and irrigator for enema insertion and rinsing. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz .
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components. Use in patients with acute inflammation of the colon, with pulmonary embolism or aortic aneurysms. Use in women with high-risk pregnancy
Adverse effects:	No adverse effects are known if the catheter is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

The catheter is intended to examine rectum, administer medicines and drainage.

Instruction for use

1. Check the type labelling of the catheter, usable life on the package and the package for damages.
2. Take the catheter out of the packaging under aseptic condition.
3. Use the catheter for the intended purpose.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	DEHP	Does not contain DEHP phthalate
Sterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	MD	Medical Device
Date	Date of Manufacture						

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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