# CHARISMA® flow

**CHARISMA® flow**

- DE Gebrauchsanweisung
- GB Instructions for use
- FR Mode d'emploi
- ES Instrucciones de uso
- IT Istruzioni per l'uso
- PT Instruções de uso
- NL Gebruiksaanwijzing
- SE Bruksanvisning
- DK Brugervejledning

- NO Bruksanvisning
- FI Käyttöohjeet
- GR Οδηγίες Χρήσης
- CZ Návod k použití
- HU Használati utasítás
- LV Lietošanas instrukcija
- LT Naudojimo instrukcija
- PL Instrukcja obsługi
- RU Инструкция по применению

**Heraeus**

## CHARISMA® flow

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## CHARISMA® flow

Charisma flow is a flowable, light curing, radiopaque hybrid composite used for adhesive, tooth-coloured anterior and posterior restorations. Due to its low viscosity Charisma flow provides excellent coverage of the tooth structure.

**Composition:** Charisma Flow is based on multifunctional methacrylate monomers (EBADMA/TEGDMA) and contains approximately 62% by weight or 38% by volume inorganic fillers such as Ba-Al-F-silicate glass and SiO<sub>2</sub>. The filler particle size is between 0,005 µm and 5 µm.

**Shades:** A1, A2, A3, A3,5, B2, B3, OA2, OA4, OB2, SB, Baseline

**Applications:**

- Fissure sealing
- Envelope fissure sealing
- Class V fillings
- Minimal invasive Class I and Class II fillings in areas not subjected to masticatory forces
- Minimal invasive Class III fillings
- Smaller contour and shade adjustment on the enamel and dentine
- Small surface repairs to direct and indirect restorations in combination with a suitable bonding agent
- Small surface depth
- Cavity lining
- Bracket retainer

**Kind of application:** Preparare la preparazione per trattamento by cleaning it with a fluoride free polishing paste. Before treatment, while the teeth are still moist, select the Charisma flow shade using the Charisma shade guide.

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4. Polymerizzare each layer according to the table below using a suitable light-curing unit (e.g. Translux® Power Blue®) for 20 or 40 sec. A dispersion layer forms during polymerisation. This layer should not be touched or removed if other composite layers are to be applied.

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**Colours:** A1, A2, A3, A3,5, B2, B3, OA2, OA4, OB2, SB, Baseline

**Indications:**

- Fissure sealing
- Envelope fissure sealing
- Filling of class V
- Minimal invasive fillings of class I and II in the zone which is not subjected to masticatory force
- Minimal invasive fillings of class III
- Small corrections of colour and shape on enamel and dentine
- Repairs of small extensive fillings of direct and indirect restorations in combination with an appropriate bonding agent
- Cleaning of mobile teeth
- Enamel sealing
- Bonding of orthodontic brackets

**Mode of use:** Prepare the restorative dentine with a fluoride-free polishing paste. Select the Charisma flow shade using the Charisma shade guide.

1. Remove the carious tissue and prepare the cavity in the maximum respect of the dentine. Wash and dry the cavity. It is recommended to use a dental air/water syringe.

2. Before applying the adhesive system (for example, iBond® Self Etch, GLUMA® Comfort Bond + Desensitizer, GLUMA® Solid Bond), it is recommended to avoid contamination with moisture, blood or saliva.

3. Apply Charisma flow in the cavity (max. 2 mm, Baseline max. 1 mm) in the cavity. During application, it is recommended that the flow is slow and uniform to prevent air bubbles from forming. Use a dental air/water syringe to clean the cavity.

## CHARIS